


Product Manual: Adapted Limit Switch



Revision: 1.0

Date: 03-Dec-2025

Table of Contents

Table of Contents.....	1
1. Introduction.....	2
1.1 About the Community Partner.....	2
1.2 Product Description.....	3
WHY?.....	3
WHO?.....	3
WHAT?.....	3
HOW?.....	3
1.3 Product Overview and Features.....	4
1.4 Bill of Materials (BOM).....	5
2. User Guide: Instructions for Use.....	6
2.1 Set Up Instructions.....	7
2.2 Maintenance.....	10
2.3 Troubleshooting.....	10
3. Maker Guide: Assembly Instructions.....	11
3.1 Assembly Guide.....	12
4. Credit and Open Hardware Licence.....	15
4.1. Open Hardware Attribution.....	15
4.2. License Statement and Source Availability.....	16
4.3 CERN Open Hardware Licence.....	16

1. Introduction

1.1 About the Community Partner

Our journey began in 1970 with the will of Mrs Shakuntala Bhatia. Under the banner of the Asian Women's Welfare Association, Mrs Bhatia and her peers, like Mrs Tambyah and Mrs Kula, identified gaps in the community to help the underserved.

The Association incorporated AWWA Ltd on 7 January 2015, a Singapore company limited by guarantee. In April of the same year, AWWA Ltd took over the operations and activities previously managed by the Association. In 2022, the Association was dissolved, and AWWA continues to address social service gaps as they arise, by being guided by our mission and strategic vision.



1.2 Product Description

This table clarifies the intended purpose and scope of the product:

WHY?	Special needs students with low mobility are often excluded from using toys and/or face lack of appropriate ways to express themselves.
WHO?	For special needs students with low mobility who are learning to understand cause-and-effect and have fun.
WHAT?	A limit switch with an additional 3.5mm male mono audio cable to trigger the power of other toys or play back audio. 3D printed casing with camera mounting hole for easy install clamp and mount on the wheelchair
HOW?	This product transforms an off-the-shelf limit switch into a durable and responsive device that triggers power on/off other toys or playback audio and an easy mounting mechanism . It allows students to use the limit switches to experience and understand the cause-and-effect.

1.3 Product Overview and Features



1. Adapted Limit Switch support push in all directions to trigger. Left, right, up or down it can be triggered.
2. Add-on soft ball, soft bicycle handle or soft toy to trigger stick to make it more interesting to students to try it out.
3. Can support different types of mounting or clamp system with standardised camera mount.
4. Connect to other electronic toys with a female 3.5mm mono audio jack cable.

<https://www.engineeringgood.org/bespoke-projects/>

1.4 Bill of Materials (BOM)

Item	Supplier/Product Link	Price	Quantity
3D printed top part	Pending for stl file		1
3D printed bottom part	https://3dprintsingapore.com/	2.77 SGD/pcs	1
Limit switch	Click Here - Link to long	22.08 SGD	1
3.5mm male mono audio jack cable	Click Here - Link too long	5.42 SGD/10pcs	1
M3 x 20mm partially threaded screw	Click Here - Link too long	6.16 SGD/3pcs	2
M3 x 6mm screw	Click Here - Link too long	1.80 SGD/pcs	1
M4 x 16mm bolt	Click Here - Link too long	1.67 SGD/5pcs	2
M4 x 30mm bolt	Click Here - Link too long	1.90 SGD/5pcs	2
M4 nut	Click Here - Link too long	1.87 SGD/20pcs	4
M5 Tee nut	Click Here - Link to long	1.41 SGD/10pcs	1
Bicycle soft handle grip	https://shop.daisosingapore.com.sg/products/4550480051574	2.18 SGD/2pcs	1
Soft ball	Click Here - Link to long	0.73 SGD/pcs	1

2025/2026

Adapted Limit Switch Product Handbook

Rev: 1.0



2. User Guide: Instructions for Use





<https://www.engineeringgood.org/ bespoke-projects/>

© 2024 Engineering Good

2.1 Set Up Instructions

Use the following table for a step-by-step guide to setting up the product:

Setting up	
Push the white or spring to the left or right to trigger it.	
Connect the 3.5mm male mono audio jack cable to toys	

Setting up

Suction cup mount. Secure the device properly by pushing the locking lever at the bottom side of the mount. It should be used on a flat surface.



Clamp mount. Secure the device properly by turning the big knob and change the angle with the small knob on the mount. It is mostly used on a wheelchair.



Setting up

Gooseneck mount. Secure the device properly by squeezing the clip part and the gooseneck can adjust the angle. It is mostly used on a wheelchair.



2.2 Maintenance

1. Check the clamp mounting if it is loose.
2. Check the 3.5mm male mono audio jack cable if it is loose or damaged.

2.3 Troubleshooting

Problem	Possible Cause	What You Can Try
The clamp or mount is loose so didn't trigger toy	Didn't clamp or mount properly	Read user guide on how those clamp and mount works
Not triggering the toy even is mounted properly	Toy issues	Check on the Toy product manual
	Limit switch connection broken	Open with the assembly instruction and try to check on it or contact EG for support.

2025/2026

Adapted Limit Switch Product Handbook

Rev: 1.0






3. Maker Guide: Assembly Instructions








<https://www.engineeringgood.org/bespoke-projects/>

© 2024 Engineering Good

3.1 Assembly Guide

Steps	Images
Remove all the 3 screws from the limit switch.	
Put the 3.5mm male mono audio jack cable through the flexible strain relief.	
Unscrew the top 2 "NO" screws until the wire from the 3.5mm male mono audio jack wire can go through it.	

Steps	Images
Screw back the top 2 "NO".	
With the 3D printed bottom part and insert the tee nut into the hole on the flat side and use a M5 bolt to screw so the tee nut will bite into the 3D printed part.	
With the 3D printer top and cover up the limit switch and screw the 3 screw and 4 bolt/nuts.	

Steps	Images
<p>Screw the different mount or clamp into the bottom tee nut.</p>	
<p>Test with a toy to see if it works.</p>	

4. Credit and Open Hardware Licence

This section provides the required legal notices and attribution for the open hardware design used to create this product.

4.1. Open Hardware Attribution

The core design for this Adapted Limit Switch is based on the Open Wobble Switch, an open-source assistive technology project.

- Original Designer: Makers Making Change (a program of Neil Squire).
- Original Copyright: Copyright (c) Neil Squire / Makers Making Change.
- Original Source Location: The original project source files are available online at:

<https://www.makersmakingchange.com/product/open-wobble-switch/01tJR000000698oYAA>



The screenshot shows the Makers Making Change website. At the top is a navigation bar with the logo, menu items (Assistive Devices, Make, Engage, Learn, Clubs That Care, More), and a search bar. The main heading is "About Makers Making Change" with the tagline "Empowering people with disabilities since 1984". Below this is a section titled "Our Mission" with text about the organization's goals and a video player showing three people working on a project. The video player has a "Watch on YouTube" button.

<https://www.engineeringgood.org/bespoke-projects/>

4.2. License Statement and Source Availability

The hardware design used in this product is licensed under the CERN Open Hardware Licence Version 2 – Weakly Reciprocal (CERN-OHL-W 2.0) or later. By distributing this product, we are obligated to make the complete design source available to you.

- Complete Source Availability: The complete design files (Source), including schematics, assembly instructions, and any modifications made by Engineering Good, are available free of charge at a permanent online location:
<https://github.com/Engineering-Good/T4G-Adapted-Limit-Switch>
- Modification Notice: This version of the Adapted Limit Switch design was modified by Engineering Good to include the 3.5mm mono audio jack functionality for triggering other toys. The full, modified source is available at the URL listed above.
- Full License Text: The complete legal text of the CERN Open Hardware Licence Version 2 – Weakly Reciprocal** follows this section.

4.3 CERN Open Hardware Licence

The full text of the license is available here: [CERN-OHL-W 2.0 Full Text](#)